



Victorian Inclusion Agency Specialist Equipment Library Catalogue June 2023

Inclusion Support Program (ISP) is funded by the Australian Government Department of Education. As part of the ISP, the Victorian Inclusion Agency is led by Community Child Care, and is delivered in partnership with Yooralla and KU Children's Services.







Inclusion Support Programme Guidelines

2016-2017 to 2018-2020

5.4. Accessing specialist equipment through the Specialist Equipment Library

The service's need for specialist equipment (<u>Section 2.1.3</u>) must be identified through the development of its Strategic Inclusion Plan, or recommended by a relevant professional, such as an occupational therapist or physiotherapist, with knowledge of the child's biomechanical functional needs. If equipment is recommended by a relevant professional, the Inclusion Agency must endorse and agree to the recommendation before an application is submitted.

A Strategic Inclusion Plan must be completed by the service in receipt of equipment through the Specialist Equipment Library within a set timeframe, as agreed by the Inclusion Agency. This means that specialist equipment may be delivered to the service before a Strategic Inclusion Plan is completed, and this consideration will form part of the loan agreement with the Specialist Equipment Library, as determined by the Inclusion Agency.

It is the responsibility of each service accessing equipment from the Specialist Equipment Library to ensure it is fitted by an appropriate professional. In some cases, an occupational therapist or physiotherapist who has knowledge of the child's biomechanical functioning and needs may be required.

Before accessing equipment from the Specialist Equipment Library, the service must sign a Loan Agreement Form with its Inclusion Agency, which sets out the Terms and Conditions associated with borrowing equipment from the Specialist Equipment Library. The process for accessing equipment from the Specialist Equipment Library, including expected timeframes, transportation arrangements, return arrangements and the Terms and Conditions associated with borrowing the equipment are determined by each Inclusion Agency. The Inclusion Agency or inclusion professional will follow-up with the service after receiving the equipment to ensure that the equipment meets the desired needs of the service (i.e., addresses the barriers identified in the Strategic Inclusion Plan) and the equipment is being used appropriately and correctly by educators within the service.

The service can access information about the Specialist Equipment Library directly from its Inclusion Agency. Contact details for the Inclusion Agency in each jurisdiction are available on the <u>department's website</u>

Version 2.1 March 2020

Content's page

1.	Activity Chairs and Seating	Page 4
2.	Floor Sitters and Positioning Systems	Page 19
3.	Standing Frames	Page 30
4.	Tables and Desks	Page 44
5.	Toileting and Change Tables	Page 49
6.	Manual Handling Aids and Appliances	Page 56
7.	Physical Activity Devices	Page 59
8.	Bead and Mirrored Stands	Page 64
9.	Switch Activated Toys Communication Devices and other Toys	Page 67
10.	Key Word Signs and Auslan Based Resources	Page 80
11.	Environment Aids and Appliances	Page 85
12.	Switches and Accessories	Page 88

ALL EQUIPMENT SHOWN IN THE CATALOGUE IS A GUIDE ONLY AND IS DEPENDENT ON AVAILABILITY AT THE TIME THE EQUIPMENT REQUEST IS ASSESS BY THE LIBRARY'S THERAPIST.

FOR PRESCRIBING THERAPISTS - EQUIPMENT SIZES AND ACCESSORIES CAN BE ACCESSED VIA THE MANUFACTURERS WEBSITES. ALL EQUIPMENT APPLICATIONS SHOULD LIST SIZES AND ACCESSORIES REQUIRED.

1. Activity Chairs and Seating

High Complexity Chairs

The Atom is the ideal seating system to give children with additional postural needs the right amount of support to develop correct seating posture and function. As children grow and develop, all support features are easily adjusted, added or removed.

The Atom by Jenx has a wide range of accessories and optional extras, meaning it can be highly customised for every child's needs.

Jenx Atom Seating System – Wooden Base or Y-Base Age Range: Approx. 6mths-5 years

Back Angle (posterior-anterior): 15 degrees-15 degrees Back Pad Height from Seat to Top: 270mm-380mm

Chest Width: 200mm-280mm Hip Width: 180mm-270mm Seat Depth: 90mm-340mm

Seat Depth with Long Leg Conversion: 315mm-565mm Lower Leg Length with Standard Footplate: 165mm-325mm Lower Leg Length with Flip-up Footplate: 115mm-360mm Footplate Angle (dorsi-plantar): 5 degrees-35 degrees



A fully adjustable, modular seating system specifically designed for young children from 4 months through to 5 years.

Jenx - Bee Seating Unit

The Bee seating system is an ideal introduction to supported functional seating.


It provides an outstanding level of versatility via fixed and height adjustable base options, along with a wide variety of head, trunk, leg and foot support options.



The giraffe chair is suitable for a child who needs extra seating support.

Jenx - Giraffe Chair It can be adjusted to accommodate a child from approximately 12 months to 5 years, and includes a range of adjustable positioning supports.

The giraffe chair can be used with a tray,or pushed up to a table without a tray (depending on table height) to enable participation in group activities.



The Jenx Junior Chair is a fully adjustable contoured modular seating system.

It has tilt in space with a hydraulic or powered base.

The system provides excellent postural support for a wide range of users.

For prescribing Therapists, please refer to the Jenx website for details on sizes and accessories available.

The advance chair is suitable for a child who needs maximum seating support. It includes a range of adjustable positioning supports. Leckey -The advance chair is height-adjustable and can be used with a tray or pushed up to a table without a tray to enable participation in group Advance Chair activities. For prescribing Therapists, please refer to the Leckey website for details on sizes and accessories available.

Leckey - Everyday Seat	Designed specifically for children or young adults with moderate needs who still require pelvic stability but do not need the same advanced level of support as those with complex needs. For prescribing Therapists, please refer to the Leckey website for details on sizes and accessories available.	
MediFab Discovery Shuttle Seat on a Launcher Base	The Launcher Base has been designed to be used with the Spex Discovery Seating System and can be raised for table height activities or lowered for floor level activities. Use the tray (optional) for feeding or play. Ultra-lightweight & compact. Tilt-in-space seating. Recline adjustment. Growth-adjustable. Wipe-down upholstery. For prescribing Therapists, please refer to the MediFab website for details on sizes and accessories available.	

Moderate Complexity Chairs

The Brookfield is a robust chair for children from the age of 2 years who require moderate levels of postural support when seated.

The Brookfield has a unique moulded back which promotes central alignment and gives additional support whilst your child is seated.

Brookfield Chair with Mobility Base Its seat height and depth can also be adjusted to suit. As an additional feature, the back can be adjusted in height.

The Brookfield comes complete with a moulded back, contoured seat, pelvic strap and height adjustable arms as standard.

A static or mobile version is also available as well as a choice of high or low back. Available as a static or mobile chair.





Rifton Compass Chair	The contoured seat and sturdy armrests provide reassuring boundaries, giving mildly involved students the support they need to feel secure and comfortable. The five available sizes accommodate everyone from small children to adults, and the adjustable legs allow fine-tuning for optimal positioning.	
Onan	For prescribing Therapists, please refer to the Rifton website for details on sizes and accessories available.	

	The Easy Seat is suitable for a child who needs a kinder chair with extra support. It has adjustable armrests, footplate, and pommel, and additional side supports are available.	
Leckey - Easy Seat	It is height-adjustable to enable a child to sit at a kinder table and participate in group activities. For prescribing Therapists, please refer to the Leckey website for details on sizes and accessories available.	
DES - Kelly Chair	The Kelly chair is suitable for a child who needs a kinder chair with extra support. It has adjustable armrests, footplate, and pommel and tray. Available is pre-school and primary sizes, each of which can be adjusted for growth. It is height-adjustable to enable a child to sit at a kinder table and participate in group activities.	

With multiple options to cater to a variety of users, the Zoomi is simple to configure. Both the seat and footrest are easy to adjust to the appropriate height and depth allowing the chair to accommodate to the child's needs.

Smirthwaite Zoomi Seat

Standard Zoomi chair comes complete with flat seat and backrest and pelvic strap as standard.

Accessories for this chair are - low back padded cover, armrest, quick release tray, 5-point harness, flat padded seat, slotted footrest and mobile skis, it offers the ideal solution for children who require a little extra postural support.

Please refer to the Smirthwaite website for sizing and accessories



The Star X Chair is suitable for a child who needs extra seating support.

It will accommodate a child from approximately 12 months to 4 years.

It includes a range of adjustable positioning supports including thoracic, hip, and foot supports.

Star X Chair

The Star X Chair can be moved up to a table (depending on table height) to enable participation in group activities.



MediFab – Neat Zeat Seating System	The Neat Zeat offers a modular indoor seating solution for activity, feeding and learning opportunities. Constructed of durable hardwood with a lacquer finish. Height adjustment grooves on the frame provide a huge growth range that can be easily adjusted with an Allen key. Three packages are available, with the addition of further accessories to script according to individual requirements. Neat Zeat is previously known as the Zeat Seat. For prescribing Therapists, please refer to the MediFab website for details on sizes and accessories available.	
Rifton Activity Chair – Alternative configurations	The Rifton Activity Chair offers a wide variety of features and accessories including a hi lo base, making it an option for most children with special needs. No tools required for adjustments, make adjustments while the child is sitting in the chair. Tilt-in-Space is a built-in feature offering 15 degrees of forward tilt and 15 degrees of backward tilt. For prescribing Therapists, please refer to the Rifton website for details on sizes and accessories available.	

Low Complexity Chairs

Leckey - **PAL** Childrens Seat An adjustable chair suitable for children with mild to moderate postural support needs. Features a wooden frame, angle adjustable backrest, adjustable seat depth, height adjustable legs, footplate, armrests, and bright upholstery.

Available in a range of sizes, please refer to the manufacturer's website for this information. A range of accessories are also available including footplate, side pads, mobile and static skis, ramped cushion and pommel.



A wooden framed chair with a padded contoured vinyl seat and backrest.

The backrest features adjustable lumbar support and recline.

The seat height and depth are adjustable.

Smirthwaite -Juni Chair

The armrests are made of wood and are height adjustable. A range of accessories are available including side pads in 25mm and 50mm, ramped cushions, harness, footplate and pommel.



Special Tomato Soft Touch -Booster Seat with Liner Back The Special Tomato Soft Touch Booster Seat with Liner Back is comfortable yet supportive seating solutions that are Latex Free, seam free, peel-resistant and impermeable to fluids. They are soft to touch for comfort, but also durable and hygienic as they are made of anti-microbial material. The Booster seats are easily cleaned with common disinfectants and so are great for use with multiple children.

Assist in promoting proper head, neck and trunk alignment

Can be attached to a standard chair with attachment straps

Can be attached to the optional Tilt wedge or Mobile tilt wedge



2. Floor Sitters, Positioning Systems, Accessories and Other:

Floor Sitters, Positioning Systems, Accessories and Other: HIGH **COMPLEXITY** The Firefly GoTo Seat by Leckey is a lightweight, portable, postural support floor seat for kids with disabilities. It's a safe, supportive multipurpose seat, which provides additional upper body postural support for children with disabilities. Leckey - Firefly Sitting securely and independently in the Firefly GoTo Seat gives GoTo Seat with children the opportunity to develop lots of skills including gross motor, Floor Sitter and fine motor, visual motor, sensory, social and communication. Pressure Cushion Pressure Cushion is also available. For prescribing Therapists, please refer to the Leckey website for details on sizes and accessories available.

A foam chair which utilises gravity to hold an individual with low postural control in a semi upright position. The chair features raised armrests to support the user laterally. The seat cover which is incontinent resistant can be removed for cleaning. Standard and mobile tilting base options are available. For prescribing Therapists, please refer to the MediFab website for details on sizes and accessories available.	HIGH COMPLEXITY		
	Gravity Chair	with low postural control in a semi upright position. The chair features raised armrests to support the user laterally. The seat cover which is incontinent resistant can be removed for cleaning. Standard and mobile tilting base options are available. For prescribing Therapists, please refer to the MediFab website for details on sizes and accessories	

HIGH COMPLEXITY		
Jenx - Corner Seat	The corner seat is suitable for a child who needs support to sit straight legged on the floor. It is available in three sizes, each of which can be adjusted for growth. The corner seat can be used with the nursery table to enable access to toys. Nursery tables are also available. For prescribing Therapists, please refer to the Jenx website for details on sizes and accessories available.	

HIGH COMPLEXITY

The Tumble Forms Floor Sitter is suitable for a child who needs extra support for floor sitting and cannot tolerate a fully upright position.

Tumble Forms Floor Sitter and Easel Table It is available in three sizes (the range accommodates children from approximately 6 months to primary school age), and the seat can be positioned on the base in varying degrees of recline.

The Tumble Forms Floor Sitter can be used with an easel table, as shown to allow access to toys and activities.



MODERATE COMPLEXITY

Scallop Seat

The Scallop is suitable for children who require a little extra help while floor sitting or sitting on a chair.

This product is also great for children with learning or behavioural disorders as it provides tactile stimulation which helps to reassure the child and enable them to relax and concentrate while using the product.



LOW COMPLEXITY Bilibo Seat	Deceptively simple and magically attractive, Bilibo's iconic shell shape and bright colours immediately appeal to children. Bilibo arouses curiosity, engages the imagination and playfully trains basic motor skills and balance. Both boys and girls find unlimited ways to play with Bilibo – using it to rock and spin in, hide under, sit on, tote with and peek through. Bilibo comes in only one size.	
LOW	Bumbo chairs can be used on any flat ground level surface. They are a lightweight, portable plastic seat suitable for babies four months and older.	
Bumbo Chair	It is shaped to support the baby when sitting. Available in only one size. The baby must have independent head control before using this product.	
	Equipment measurements: - H 20CM x W 13CM x L 13CM	

LOW COMPLEXITY		
	Lightweight and portable	
Karibu Seat with Tray	To support unstable sitting posture. The baby must have independent head control before using this product.	
, , ,	Easy to care	
	Stains can be easily wiped off with a wet towel	
	Size: L44.5 x W34.7 x H23cm	South
LOW COMPLEXITY	Four independently adjustable non-slip feet let you sit on any chair.	
	This seat has a lightweight body with a carry handle.	
Mamas & Pappas	Great as they grow - soft and comfortable padded insert for use from 6 months to 3 years.	
Booster Seat	Extra safe - two chair harnesses' make sure they stay put.	50
with Tray	Age Suitability: 6 months or when a child can sit unaided - 3 Years Approx.	988

LOW		
Move'n' Sit Inflatable Cushion	An inflatable wedge which serves as a seat wedge and encourages an upright posture. The upper surface is bumpy and provides tactile feedback while undertaking other activities.	

Positioning System

This lightweight bag contains an entire activity centre.

Inside you will find rolls, wedges, and supports (designed by clinical experts) that stick together in dozens of different ways to let you and your child practice the early positions which form the building blocks for later abilities.

Leckey Firefly Playpak

Equipment measurements in bag: - H 20CM x W 45CM x L 45CM



Positioning System

Allows patient to be positioned in prone - good positioning for eliminating the effect of gravity present when individual is upright.

Wedge with Detachable Headrest Allows patient to be positioned to stretch flexed hips and knees. Scoop shape decreases chance that child could roll off wedge.

Headrest is useful for child with poor head and neck extension.

Equipment measurements: - H 30 CM (at high end of wedge) x W 43CM x L 150CM



Accessories and Other:		
Jenx - Nursery Tables	The Jenx Nursery Table is available in two sizes to work with the Jenx Corner Seats. It's a sturdy table with a large, clear surface area that is big enough for a child to play or work at most activities with ease. The table (and its surface) are ideal for drawing and writing activities. For prescribing Therapists, please refer to the Jenx website for details on sizes and accessories available.	

3. Standing Frames

Multi-Standers, Sit-to-Stand Standing Frames, Supine Standers, Prone Standers and Upright Standers

EasyStand – Supine / Full Standing Bantam Stander / Chair The Bantam line of standers are the only sit to stand standers that can combine the option of supine positioning to create infinite positioning possibilities from 90-90 sitting, fully vertical standing, full supine lying and anywhere in between.

No other standers can provide the amount of positioning that is possible with the EasyStand Bantam.

For prescribing Therapists, please refer to the The Altimate Medical Industries (AMI) website for details on sizes and accessories available.



Jenx -Multistander – Supine / Prone / Upright Standing System This standing frame is suitable for children between the ages of nine months and six years old.

It can be configured to support supine, prone, and upright standing with two types of headrests.

The Multistander also allows for a tray to be attached in the supine, prone and the upright position.



A versatile standing system offering Prone or Supine standing that accommodates abducted or neutral leg positioning for children between 1 – 9 years. With an extensive range of accessories, Standz can accommodate children with a range of postural requirements.

Abduction
Stander

Age Range: 1 to 9yrs Angle 0 - 90°

Age Range: 1 to 9yrs Angle 0 - 90°

For prescribing Therapists, please refer to the Jenx website for details on sizes and accessories available.



Three-in-one stander, offering supine, prone and upright positioning. Ergonomically designed chest pad replicates the shape of the rib cage. Fixed hip and knee contractures can be accommodated via the Leckey Mygo mechanically linked movement of the swinging footplate. Stander For prescribing Therapists, please refer to the Leckey website for details on sizes and accessories available.

The Squiggles Stander is a bright colourful three-in-one stander for younger children.

Leckey -Squiggles Standing Frame – Prone / Upright / Supine This versatile, adjustable stander offers prone, upright, and supine standing in one product, allowing the therapist to offer standing therapy to address a child's needs and support inclusion.



The DES Supine Standing Frame is a supine (reclined)/upright standing frame suitable for a child who needs support to stand and may not be able to tolerate a fully upright position and includes a range of adjustable positioning supports.

DES - Supine Standing Frame

It can be adjusted to enable a standing position from horizontal to fully upright.

The DES supine standing frame has an activity tray to allow access to toys.

Equipment measurements: - H 106CM x W 40CM x L 70CM



The Jenx Supine Stander is a robust stander with comfortable, adjustable accessories to provide safe standing to suit children with a wide variety of needs. While supported and upright, children can participate in a range of activities.

Jenx - Supine Stander Available in two sizes to suit children from 9 months onwards.

A robust stander with a range of accessories including trays to suit a variety of needs.



The Humphrey is a prone (forward-tilted) upright standing frame suitable for a child who needs support to stand and may not be able to tolerate a fully upright position.

Jenx – Prone -Humphrey Standing Frame It can be adjusted to accommodate a child from approximately 3 to 7 years and includes a range of adjustable positioning supports. It can be adjusted to enable a standing position from horizontal to fully upright.

The Humphrey standing frame has an activity tray that can be used either with an insert to create a flat play surface, or with a bowl for sensory play.



The Monkey is a prone (forward-tilted) upright standing frame suitable for a child who needs support to stand and may not be able to tolerate a fully upright position.

It can be adjusted to accommodate a child from approximately 12 months to 5 years and includes a range of adjustable positioning supports.

Jenx – Prone -Monkey Standing Frame It can be adjusted to enable a standing position from almost horizontal to fully upright.

The Monkey Standing Frame has an activity tray that can be used either with an insert to create a flat play surface, or with a bowl to allow sensory play.



The Jenx prone standing frame is a prone (forward-tilted) upright standing frame suitable for a child who needs support to stand and may not be able to tolerate a fully upright position.

It can be adjusted to enable a standing position from upright to 40 degrees.

Jenx - Prone -Standing Frame

The Jenx prone standing frame has an activity tray and bowl to allow access to toys.



The Jenx Rabbit Paediatric Mobile Stander is a secure, reliable, and easy-to adjust standing aid that supports special needs children in the vertical and prone positions.

Jenx – Prone -Rabbit Stander The tray comes with a bowl insert.

Equipment Measurements: - H 60CM x W 38CM x L 55CM



This prone standing frame features multiple adjustments on all components to grow with the child.

It consists of a wooden frame, padded vinyl supports, and quick release straps with buckles.

DES – Prone -Standing Frame

Frame incorporates a manually operated tilt adjustment with a positive locking mechanism and a solid back support to hold child in position.

Mobile frame folds flat for storage.

Equipment measurements: - H 96CM x W 55CM x L 73CM



The Dynamic Standing Frame is suitable for a child who needs support to stand and is designed to allow the child some movement.

Healthtec
Dynamic /
Chattanooga –
Upright
Standing Frame

It is available in the small and medium size and can be adjusted in height. Another option is the Standing Frame can be attached to a platform with wheels for easier mobility.

The Dynamic Standing Frame can be used with a Taddy height adjustable cut out table, or it can be pushed up to an easel, water trough or high table to allow access to activities with other children.





Jenx – Upright -Penguin Stander The Jenx Penguin Stander has been designed to provide excellent support for children in upright standing, and the range and support of adjustment options available creates a product that is suitable for children aged 18 months to 4 years.



4. Tables and Desks

Tables and Desks

DES - Height Adjustable Easel Table The height adjustable easel table provides a slanted surface for activities such as drawing and reading.

It is angle-adjustable, bringing activities into the child's line of vision, improving posture and encouraging fine motor development.

Equipment measurements when collapsed to its lowest setting -

H 26CM x W 62CM x L 42CM



Height Adjustable Table with Fold Down Pin Legs	Height adjustable with fold down legs to provide flexibility with ensuring the table is at the correct height for the child. For prescribing Therapists, please refer to the CAP Furniture website for details on sizes and accessories available.	
Height Adjustable Table (with wheelchair cut out)	Height adjustable table with left or right hand cut out for power wheelchair control. For prescribing Therapists, please refer to the CAP Furniture website for details on sizes and accessories available.	

Height Adjustable The Library stocks a variety of table sizes.

These tables are engineered to be height adjustable with laminated tops and steel frames.

Square Tables are useful in a variety of settings and the most standard inclusion in OSHC school programs and Childcare Centres.

The Library stocks a variety of table sizes.

For prescribing Therapists, please refer to the CAP Furniture website for details on sizes and accessories available.

Leckey PAL Table Non Adustable - PINK	Pink Non-adjustable table. Equipment measurements: - H 51CM x W 50CM x L 130CM	
Taddy Table	Height adjustable legs with table cut out provides flexibility when used in conjunction with standing frames or sitting systems. Equipment measurements: - H 58CM W 50CM L 50CM	

TAD –
Adjustable Table is suitable for young children who need support to maintain an upright and stable posture.

Table width 600mm Table depth 500mm Cut out width 220mm Cut out depth 95mm

Table height from ground* 230-455mm - 230-455mm

*Lowest possible heights are achieved with adjustable legs removed.

For prescribing Therapists, please refer to the TAD Adaptable Solutions website for details on sizes and accessories available.

5. Toileting and Change Tables

Toileting		
Electric Change Table	The electric change table is suitable for a larger child who is able to do a standing transfer to and from the change table, reducing the need for lifting. It has an electric height-adjustment mechanism to enable it to be lowered for transfers and raised to a suitable working height for nappy changes. Cot sides or rails provide extra security. The electric change table is fitted with castors to enable ease of movement for storage. Various sizes stocked in Library.	
Flamingo Toilet Chair	The R82 Flamingo Toilet chair provides extra toileting support. It is available in several sizes to accommodate children from approximately 3 years to adolescent. It has a range of adjustable positioning supports and can tilt in space. It can be used as a free-standing commode or positioned over the toilet. Please refer to the manufacturer's website for sizing and accessories	

Platform Step and Rails

This multi-purpose step up platform has support handrails on both sides to give 240mm extra height.

It is designed to step up to toilets, washbasins, tables, and other areas where a sturdy, stable step is required. It has a slip resistant surface.

Equipment measurements: - H 55CM x W 51CM x L 75CM

Height adjustable platform with slip resistant mat and adjustable handrails.

Contemporary design with comfortable, smooth handles

Nuvo Petite Toilet Platform Lightweight, powder coated steel frame makes it sturdy yet easy to move

Non-slip rubber ferrules on the feet ensure additional safety



Potty Infant Chair - Adjustable

Contoured fibreglass mounted seat, mounted on zinc treated, epoxy coated steel frame with optional chest strap available.

Fits most toilets when fitted with adjustable legs.

The Blue Wave Toileting System is suitable for a child who needs extra support for toileting.

Rifton Blue Wave Toileting System

(Small and Large Chair Systems) It is available in two sizes to accommodate children from approximately 3 years to adolescent, and has a range of adjustable positioning supports. It can be used with or without a mobile frame.

With the mobile frame it can be used as a free-standing commode with a potty attachment, or as an over-toilet commode without the potty attachment.

Without the mobile frame the seat can be mounted on the toilet.



Supplied with a harness for full trunk support.

Smaller model is intended for children up to 10 years.

Toilet Support -High Backed

Reduced seat depth (RSD): sits 4" forward to lessen seat depth when attached to U-Shape Mounting Base.

Dimensions of supports: height adjustable 19-24 inches; width at top 8" and width at bottom 12".



Toilet Support -Low Back

Designed to fit any school toilet, the toilet supports strong and secure frame allows it to be easily adjusted for maximum user control and comfort, and also allows the low back toilet support to be quickly removed when not in use for easy storage.



6. Manual Handling Aids and Appliances

Manual Handling Aids and Appliances

The mobile hoist is suitable for situations where manual lifting of a child is not viable, usually for a larger child who is not walking and is unable to do a standing transfer.

Hoists – Battery Operated -Rechargeable It is useful for lifts across short distances, usually to and from specific positions, for example from chair to change table or chair to floor.

Slings are available in a variety of styles depending on the amount of support required.

The mobile hoist is powered by a rechargeable battery and is operated with a hand-held control.



Hoisting Slings	Manufactured under the highest quality control standards. The extremely durable slings come in a range of fabrics and sizes to suit the individual child when used in conjunction with an electric powered hoist to assist educators and staff with toileting needs. Various sizes stocked in Library. Please note – prescribing Therapists need to provide the brand and size of sling to service.	
Wheelie Stool	The wheelie stool is a mobile, height-adjustable stool designed to enable staff to work with children without bending or twisting. It has a wheeled base, an upholstered seat which can swivel a full 360 degrees and gas-assisted height adjustment. It enables staff to work at the level of the child for feeding, seated activities, and supported standing and walking.	

Wooden Stool

Height adjustable with screws with a vinyl padded seat and a wheeled base for educational staff.

Equipment measurements: - H 36CM x W 40CM x L 14CM

7. Physical Activity Devices

Bicycles, Full Support Swings, PMP Gross Motor Steps and Ramps

The Supportive Tricycle is suitable for a child who is unable to use a regular tricycle.

Supportive Tricycle

It can be adjusted to accommodate a child from approximately 3 to 6 years and includes a range of adjustable positioning supports.

The Library stocks a variety of sizes.



Full Support Swings	The Full Support Swing seat is suitable for a child who needs a swing seat with extra support. It is available in three sizes (the range accommodates children from approximately 2 years to adult) and provides trunk and head support. The full support swing seat will fit onto standard "pigtail" swing hooks and can be adjusted in height and angle of recline.	
------------------------	---	--

JennSwing semi-reclined design offers upper and lower body support with added leg rest.

The arm rests provides lateral support. Includes easy adjustable 5 point safety harness that secures the child by both shoulders and there waist.

Full Support Swing Seat -JennSwing Adaptive Swing The seat is 14" wide x 11" deep.

The seat back is 24" high and knee to foot length is 14" long. It is 45" overall and has a 3/16" thickness.

This special needs swing seat is for children weighing up to 56.7kg.

There is a fully adjustable safety harness that securely locks the child into place.

Supplied with 7ft chain length.

Attaches to most swing sets quickly and easily with S-hooks.

Please note measurements are in inches.



Glentham - Vinyl Wedges and Shapes	The Library has a variety of sturdy and waterproof equipment that can promote inclusion and play activity. Types of equipment are wedges, various sized cylinders, and different shaped animals including Loch Ness and Elephant, soft-sided See Saw and double sided – 2 steps.	
Fixed Wooden Step	Equipment measurement: - Fixed height L 30CM x W 37CM x H 30CM	

Platform Steps	Platform to help assist children with posture when sitting and accessing benches. The Library stocks a variety of sizes.	
----------------	---	--

8. Bead and Mirrored Stands

Bead and mirrored stands

Provides stimulating feedback! Specially-designed for the individual who requires exceptional sensory input, this 13cm x 20cm bead chain and unbreakable plastic mirror plays music, small red light is activated and the unit vibrates when touched.

Lighted Sensory Plate / Switch

The Sensory Plate can also be used as an external switch without vibration, music and/or lights when connected to a toy or device.

Size: 30cm x 21cm x 7.5cm. Requires 2 C Batteries

For prescribing Therapists, please refer to the Enabling Devices website for details.



Peek-a-Boo Mirror Switch adapted sensory toy	Encourages exploration! Mirror provides tactile and visual reinforcement. Playing with the beads will activate lights, vibration and music, while the soothing chains provide unique tactile stimulation. Vibration sensitivity is adjustable. Can be used as a switch to activate devices. Can be used on the floor or a table. Lights, vibration and music can be turned on or off by hidden switches. Requires 6 size AA Batteries. For prescribing Therapists, please refer to the Enabling Devices website for details.	
Bead Chain - Music Switch	Hours of tactile and auditory stimulation! This innovative Bead Chain is great for developing the concept of cause and effect. Movement of the beads will turn on a built-in music box. As long as the beads are making contact with the metal bar, the music stays on. Can also be used as an external capability switch without music when connected to any toy or device. Size: 32cm H x 37cmL x 27cmW Requires 2 AA Batteries For prescribing Therapists, please refer to the Enabling Devices website for details.	

'Razzle Dazzle' Bead Chain OR Multi-Coloured Bead Chain	The Razzle Dazzle Bead Chain Curtain provides hours of somatosensory stimulation. Who can resist running their fingers through our silver or multicoloured bead chain curtains? Mounted on a sturdy base, rubber suction cups hold the bead chain in place on almost any surface. Delivers a wonderful tactile experience for all! Size: Size: 45cm L x 20cmW x 32cmH For prescribing Therapists, please refer to the Enabling Devices website for details.	
Somatosensory Bead Chain with Mirror	What's behind beaded curtain? Surprise! It's an unbreakable plastic mirror that reflects the face of the user! A terrific toy that provides visual and tactile stimulation and body identification! No batteries required For prescribing Therapists, please refer to the Enabling Devices website for details.	

9. Switch Activated Toys, Communication devices and other Toys

Switch Activated Communication Devices and Toys Specially designed for children with limited mobility. This small-scale activity centre is designed to provide sensory stimulation and cognitive activity. Great for the user whose physical mobility may be limited due to positioning needs. Features a bead curtain, a dome with brightly-coloured seguins which spin when the button is pressed, a textured plate which Compact **Activity Centre** vibrates and lights up when pressed, and a pull ball that turns on a Sensory Toy music box. Requires 2 size C batteries to operate. For prescribing Therapists, please refer to the Enabling Devices website for details.

There are 5 different switches built into the brightly coloured activity centre. Each has a different function that has been designed to stimulate sensory awareness, hand-eye coordination and aid in teaching cause and effect:

- The bright yellow plate activates a built-in AM/FM radio.
- The large red plate vibrates when pressed.
- Moving the wooden roller sounds a buzzer.
- A pull ball turns on a music box.
- The large orange plate turns on a light.

Size: 96.5 x 71 x 10 cm

Requires: 2 size C batteries

For prescribing Therapists, please refer to the Enabling Devices website for details.



Five Function Activity Centre

Sensory Toy

Music Blocks - Musical Toy	Each colourful block plays a different phrase of music and each side of each block plays a variation of that phrase. You can just touch or mix and match and rearrange the blocks to create their very own musical compositions. Mozart's Night Music. High quality CD-like digital sound. A light comes on when a block is touched. The red button at the front can also be pressed for a musical note. Requires batteries 4 AA. For prescribing Therapists, please refer to the Small World Toys and NEUROSMITH website for details.	
Activity Dome – Sensory Toy	Dome has buttons that when pressed will vibrate and a cord that you can pull which vibrates the dome. For prescribing Therapists, please refer to the Enabling Devices website for details.	

Switch Toy - Ariel Mermaid	A switch adapted toy that provides cause and effect when switch is activated by the child. Requires batteries and an additional switch for operation.	
Switch Toy – Bongo Bright Beats 2 in 1 Musical Drum Roll	Musical Drum Roll is an infant drum for sitting or crawling babies and toddlers with music, sound effects, and lights. It features two modes of music and sound - Drum Beat Mode and Jam-Along Mode. Sensory toy also has clacker beads, and a colourful roller ball. Press your capability switch to hear up beat sounds, music & lights. Requires batteries and an additional switch for operation.	

Switch Toy Buzz Lightyear

A switch adapted toy that provides cause and effect when switch is activated by the child.

Requires batteries and an additional switch for operation.

Switch **Toy** -Motorized Glitter Roll Music Box

The Motorized Glitter Roll Music Box works with an external capability switch – please note the Library has two glitter rolls, this one pictured requires a switch the other glitter roll on page 84 is a switch.

When the switch is activated, the potpourri of sparkling shapes will tumble as the sweet-sounding music plays.

Terrific for increasing auditory and visual stimulation.

Size: 35cmL x 20cmW x12cmH.

Requires Capability Switch and 2 C Batteries which are not included.

For prescribing Therapists, please refer to the Enabling Devices website for details.



He flaps his ears and claps along to "If You're Happy and You Know it, Clap your Hands".

Switch Toy – My Little Puppy Requires 3 'AA' Batteries

For prescribing Therapists, please refer to the Cuddle Barn website for details.



The remote control has been switch adapted. Allowing drivers to control the vehicle with two easy-to-use directional buttons. It helps the user to grow their hand-eye coordination skill and imagination skill.

Switch Toy Small Ferrari
Remote Control
Car

Attach your switches to drive forward press your other switch and it turns in reverse

Requires: Car: 2 x AAA 1.5v

Controller: 3 x AAA 1.5v



Switch Toy – Nemo	A switch adapted toy that provides cause and effect when switch is activated by the child. Requires batteries and an additional switch for operation.	THE MC
Switch Toy - Pudgy Piglet	Hold down your switch and watch pudgy piglet walk and make animal sounds. Cute and cuddly, these toys will provide sensory feedback and promote play with other children. Switch is not included. Requires batteries for operation. For prescribing Therapists, please refer to the Cuddle Barn website for details.	

You can use this wonderful musical toy in four different ways: activate the attached yellow switch and the colourful bells twirl around, sounding the musical scale while the multi-coloured LEDs blink, or use a capability switch to activate the toy.

Ring Around Bells is excellent for music therapy and music time. The bells are precision tuned and provide exceptionally clear sounds.

Switch Toy -Ring Around Bells

Size: 30cmL x 35cmW x 25cmH.

Requires 4 size C Batteries.

For prescribing Therapists, please refer to the Enabling Devices website for details.



Switch Toy – Robot	A switch adapted toy that provides cause and effect when switch is activated by the child. Requires batteries and an additional switch for operation.	
Switch Toy - Textured roller	This captivating device provides a very special somatosensory experience for children. Bumps on the switch offer tactile stimulation and help hand-eye coordination. Roller switch can be used as a switch-only when you spin the drum and the unit is activated - or as a switch-toy when you spin the drum to activate lights vibrations and music. Requires 2 size AA batteries.	

Switch Toy - Thomas the Tank Engine with Remote Control	Follow Me Thomas's controller has been switch adapted for him to chug along in a large 2m circle, with sound and music to capture the child's attention. Requires 2 size AA batteries.	THOMAS GRIENDS Ways to Play! Thomas Hogically Polises the Light! Light the way and Thomas follows!
Switch Toy - Wheelie the Soft Bus	Wheelie the soft bus is switch adapted and moves gently side to side with his mouth opening and eyes lighting up singing "Wheels on the bus" Song plays for 1 minute & 24 seconds; you can press your switch to turn off midway through the tune. 23cm long x 15cm wide Requires 3 size AAA batteries For prescribing Therapists, please refer to the Enabling Devices website for details.	

Various
Weighted Toys

Weighted stuffed animals are the ultimate sensory toy.

The varied shapes, sizes and textures of our furry friends are not only appealing but they allow for weight to be distributed where it's needed most.

Children of all ages and sizes to benefit from the calming effect of deep sensory, proprioceptive, and tactile sensory input.

10. Key Word Sign and Auslan Based Resources

Kits, Books, Posters and Dictionaries

Everyday vocabulary in easy to remember signing pairs. For example - Girl / Boy / Please / Thank You / Hot / Cold

Auslan – Set of Flash Cards for Beginners These signs have been selected to provide high reward-benefit for young children when they utilise the signs, making them easily adoptable into daily vocabulary.

For prescribing Therapists, please refer to the Bilby Publishing website for details.



Auslan Baby Sign Book - ITSY BITSY SPIDER	This delightfully illustrated board book will help your child or children learn Auslan through song. The book was produced in the UK using BSL (British Sign Language). This is an Australian modified edition to include three insert stickers inside the book for the 1 sign that are different from Auslan (Australian Sign Language). For prescribing Therapists, please refer to the Bilby Publishing website for details.	Signal Straight Spider Singalord Spider Spider
	website for details.	child's Playo Annie Kubler
Auslan Baby Sign Book – MY FIRST SIGNS	This beautiful picture dictionary board book depicts over 50 commonly adopted baby/toddler signs. Each sign is clearly depicted, with descriptions of how to do the sign and colourful illustrations of the object being signed. Written in BSL, this edition includes insert stickers for Auslan variations. For prescribing Therapists, please refer to the Bilby Publishing website for details.	MILK SIGNS SILEEP Annie Kubler CAT Child's Playe CAR

Auslan Baby Sign Book – TEDDY BEAR, TEDDY BEAR!	This delightfully illustrated board book will help your child or children learn Auslan through song. The book was produced in the UK using BSL (British Sign Language). This is an Australian modified edition to include three insert stickers inside the book for the 1 sign that are different from Auslan (Australian Sign Language). For prescribing Therapists, please refer to the Bilby Publishing website for details.	Sign Teddy Bear, Teddy Bear, Teddy Annie Kubler
Auslan Books - Various Posters and Flash Cards	By teaching simple sign language to children from as young as seven months, we can help them to convey their emotions and their needs, as well as giving them a valuable head start in their social and intellectual development. For prescribing Therapists, please refer to the Bilby Publishing website for details.	Mayor Illrithday Company Managery Company Introduce of the second of t

The soft cover book is 29 \times 29 cm and includes holes in each page that gradually get bigger to reveal the animal inside her.

Auslan Puppet Set & DVD – Old Woman who Swallowed a Fly The doll has a large mouth and can actually swallow each of the bean filled animals as you go through the story.

Great to help children visualise what is happening!

For prescribing Therapists, please refer to the Bilby Publishing website for details.



The back of the song card includes detailed descriptions of how to sign all the insects, plus larger versions of the sign illustrations

The set also includes a large Green frog full hand puppet with wide opening month; and 6 insect finger puppets:

Furry Bee (9cm fury body with clear wings)

Beautiful Butterfly (9.5 cm wingspan, silk like body)

Auslan Puppet Set – Five Speckled Frogs with Song Sheet & 7 Finger Puppets

Hairy green and Yellow Caterpillar (9.5 cm long)

Dainty Ladybug (8cm wide)

The Giant Snail (10 cm tall)

The Hairy spider (8cm long with a big smile)

All finger puppets have elasticised entrances so that they stay firm on child and adult fingers.

For prescribing Therapists, please refer to the Bilby Publishing website for details.



11. Environmental Aids and Appliances

Grab Bars, S	Suction Rails and Ramps	
Stabilising Horizontal Grab Bars and Vertical Anchor Grab Bars	The Stabilising Grab Bar encourages the child to bring the arms forward to stabilise, improving posture and increasing access to activities. It consists of a horizontal or vertical bar that attaches to the tabletop with suction cups and can be used on the table or tray provided with a chair or standing frame.	
Suction Grab Rail	Suction grab rails are lightweight and portable, featuring a durable commercial grade suction cup making them the strongest in the market. Easy to install and removes without any tools.	

Depac Ramps	Portable folding ramps.	
Easy Access Portable Ramps	Easy access ramps are a lightweight, yet portable ramp designed for doorways, sliding glass doors, and raised landings. One metre in length.	

Threshold Ramp - Small and Large Sizes	A small threshold ramp designed for use in doorways, sliding glass doors and raised landings. The Library has a variety of sizes stocked.	
Steps with Rails	Designed for assisted child access to change tables. Construction is lightweight and the handrails are adjustable in height. Equipment measurements: - H 73CM W 52CM L 70CM	

12. Switches and other Accessories

Accessorie	es	
Big Keys Computer Keyboard	Big keys keyboards for computer access provide users large, high-visibility keys that are easy to see. In addition, for users who lack fine motor movements, the larger keys provide a larger strike area that is easier to access compared to standard keyboards.	1 6 2 7 3 5 4 7 6 5 6 8 7 8 9 0 Cox Cox
Campus SX and Microphone Boom with Case	Campus SX is a dedicated Phonak FM transmitter for educational applications when used in combination with Micro Boom. Professionals, teachers, educators and children benefit from maximum flexibility and user comfort.	PHONAK A V

Computer Mouse Assistance -Switch Adapted Version The BIGtrack Trackball Switch Adapted has a 3 inch trackball which makes it the largest trackball available. The large ball requires less fine motor control than a standard trackball and it is ruggedly built. It has a left and a right mouse click button located behind the trackball to avoid unwanted mouse clicks. The BIGtrack trackball also allows for 2nd Mouse Connection both the BIGtrack and a second mouse can be connected allowing simultaneous participation. The bright large 3" trackball helps to make mouse movements easy and accurate. While the oversized buttons make mouse clicks simple for those with poor motor control. For even more access this switch adapted version comes with two standard 1/8" plugs for two switches. Its left and right mouse clicks have been adapted so that you can use a switch for activation. This large trackball is also perfect for kids of all ages, especially children in preschool. The large trackball is easy on little hands and is simple to use.

The BIGtrack Switch Adapted is easy to install and works on any Windows or Macintosh computer that has a USB or PS/2 port.



The Phonak's DigiMaster 5000 innovative loudspeaker array features no less than 12 individual high-quality loudspeakers, all housed in a robust aluminium frame.

Phonak
Digimaster
5000 kit
including
Speaker with
stand and
Pass-around
Mic

The DigiMaster 5000 array emits cylindrical sound waves at the most important speech frequencies, ensuring that the required sounds are carried effectively across a large area, while generating very little speech-hindering 'late reverberation' or echo. The DigiMaster 5000 is supplied with a sturdy floor stand

The Pass-around Mic can be handed from one child to another, allowing the child who is wearing Roger Receivers to hear not only the teacher/educator but also other children's comments.



Record any series of messages directly into the Little Step- by-Step Communicator and press its activation surface for the first message. Press it again and Step-by-Step automatically plays the next message in a sequence.

Step-by-Step communicators allow children with even the most limited communication abilities to take turns in "an essential conversation."

Little Step by Step Communicator

Each message in the sequence can be as long or short as needed.

The total number of messages is only limited by the 2 minute recording time.

Requires batteries to operate equipment.

For prescribing Therapists, please refer to the Ablenet Inc. website for details on sizes and accessories available.



Switch – Big Mac Communicator	Record any single message directly into the Big Mac and press its large activation surface for up to 2 minutes of playback. Students with visual disabilities and those who require a larger activation surface will have no problem connecting with BIG Mac communicator's large, 13cm (5-inch) activation surface as they connect with teachers, educators, children, and the world around them. Requires batteries to operate switch. For prescribing Therapists, please refer to the Ablenet Inc. website for details on sizes and accessories available.	
Switch – Big Red	Large colourful surface makes this switch ideal for children with visual, cognitive, or physical requirements. For prescribing Therapists, please refer to the Ablenet Inc. website for details on sizes and accessories available.	

Switch - iSwitch	A bluetooth switch for wireless switch access to iDevices iOS 7 – iOS 11 Switch Control Works with iPad (all models including Pro / Air / Air 2 / 2-6, iPad Mini 1-3), iPod Touch 3-5, iPhone 3GS, 4 Bluetooth wireless interface - 75mm diameter) switch For prescribing Therapists, please refer to the Spectronics website for details on sizes and accessories available.	
Switch – Jelly Bean	Low profile circular switch on mounting plate for wheelchair tray or table top. Responsive across entire surface. Range of colours available. Switch can be mounted on small triangular plate, or two or four-place rectangular plate. Twist model enables switch tops to be interchanged with one of a range of colours. For prescribing Therapists, please refer to the Ablenet Inc. website for details on sizes and accessories available.	

Switch - Smiley face	Switches enable children with severe disabilities to experience cause and effect, to play and have fun, and to interact with their environment. For prescribing Therapists, please refer to the Ablenet Inc. website for details on sizes and accessories available.	
Switch – Glitter Roll Music Box	This innovative switch provides hours of multi-sensory stimulation with its roller action and built-in music box. The clear, cylindrical roller contains a potpourri of sparkling shapes. Roll it to activate the music base. Can also be used as an external capability switch without music when connected to a device - lead supplied Size: 35cmL x 17cmW x 12cmH. Requires 2 AA Batteries. Weight: 800grms. *** This device cannot be activated by a switch, it is manually rotated and when you use the switch lead it will activate another switch adapted toy ***	